

Product Data Sheet

Recombinant Human CCDC62 Protein

Catalog #	Source	Reactivity	Applications
CST1311	E. coli	H	
Description	Recombinant Human CCDC62 Protein is produced by E. coli expression system and the target gene encoding Leu131-Lys380 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CCDC62		
Alternative Names	Coiled-coil domain-containing protein 62; Protein TSP-NY; Protein aaa		
Entrez Gene	84660		
SwissProt	Q6P9F0		
Purity			
Chemical Structure	LEARNETLSN TLVELSAQVG QLQAREQALT TMIKDKDI IEAVNHIADC SGKFKMLEHA LRDAKMAETC IVKEKQDYKQ KLKALKIEVN KLKEDLNEKT TENNEQREEI IRLKQEKSC HDELLFTVER EKRKDELLNI AKSKQERTNS ELHNLRLQIYV KQQSDLQFLN FNVENSQELI QMYDSKMEES KALDSSRDMC LSDLENNHPK VDIKREKNQK SLFKDQKFEA ML		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com