

Product Data Sheet

Recombinant Human ANKRD16 Protein

Catalog #	Source	Reactivity	Applications
CST1312	E. coli	H	
Description	Recombinant Human ANKRD16 Protein is produced by E. coli expression system and the target gene encoding Arg201-Asp350 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ANKRD16		
Alternative Names	Ankyrin repeat domain-containing protein 16		
Entrez Gene	54522		
SwissProt	Q6P6B7		
Purity			
Chemical Structure	RDNCGVTALM DAIQCGHIDV ARLLLDEHGA CLSAEDSLGA QALHRAAVTG QDEAIRFLVS ELGVDVDVRA TSTHLTALHY AAKEGHTSTI QTLLSLGADI NSKDEKNRSA LHLACAGQHL ACAKFLLQSG LKDSEDITGT LAQQLPRRAD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com