

Product Data Sheet

Recombinant Human RNF220 Protein

Catalog #	Source	Reactivity	Applications
CST1315	E. coli	H	
Description	Recombinant Human RNF220 Protein is produced by E. coli expression system and the target gene encoding Ala291-Ala430 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RNF220		
Alternative Names	C1orf164; E3 ubiquitin-protein ligase RNF220; RING finger protein 220		
Entrez Gene	55182		
SwissProt	Q5VTB9		
Purity			
Chemical Structure	ATDDLHHS DR YQTF LRVRAN RQTRLNARIG KMKRRKQDEG QREGSCMAED DAVDIEHENN NRFE EYEWCG QKRIRATTLL EGGFRGSGFI MCSGKENPDS DADLDVDGDD TLEYGKPQYT EADVIPCTGE EPGEAKEREA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com