

Product Data Sheet

Recombinant Human UBXN11 Protein

Catalog #	Source	Reactivity	Applications
CST1316	E. coli	H	
Description	Recombinant Human UBXN11 Protein is produced by E. coli expression system and the target gene encoding Met21-Ser170 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	UBXN11		
Alternative Names	SOC; UBXD5; UBX domain-containing protein 11; Colorectal tumor-associated antigen COA-1; Socius; UBX domain-containing protein 5		
Entrez Gene	91544		
SwissProt	Q5T124		
Purity			
Chemical Structure	MNPGRRGIRI YGDEDEVDM L SDGCGSEEKI SVPSCYGGIG APVSRQVPAS HDSELMAFMT RKLWDLEQQV KAQTDEILSK DQKIAALEDL VQTLRHPAE ATLQRQEELE TMCVQLQRQV REMERFLSDY GLQWVGEPMD QEDSESKTVS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com