

Product Data Sheet

Recombinant Human CFAP70 Protein

Catalog #	Source	Reactivity	Applications
CST1317	E. coli	H	
Description	Recombinant Human CFAP70 Protein is produced by E. coli expression system and the target gene encoding Asn751-Leu880 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TTC18		
Alternative Names	Tetratricopeptide repeat protein 18; TPR repeat protein 18		
Entrez Gene	118491		
SwissProt	Q5T0N1		
Purity			
Chemical Structure	NNDIRMEMAF HEASKQLQAR MLQAQVTKQK STGVEDTEER GKRESSLGPW GITNGSATAI KVEAPAGPGA ALSILDKFLE ESSKLQSDSQ EPILTTQTWD PSISQKPSNT FIKEIPTKKE ASKCQDSSAL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com