

Product Data Sheet

Recombinant Human RAB41 Protein

Catalog #	Source	Reactivity	Applications
CST1318	E. coli	H	
Description	Recombinant Human RAB41 Protein is produced by E. coli expression system and the target gene encoding Met11-Leu210 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RAB41		
Alternative Names	Ras-related protein Rab-41		
Entrez Gene	347517		
SwissProt	Q5JT25		
Purity			
Chemical Structure	MEAGGFGLEA AERTEYQSLC KSKLLFLGEQ SVGKTSIISR FMYNSTFGCAC QATVGIDFLS KTMYLEDDQIV QLQLWDTAGQ ERFHSLIPSY IRDSTIAVVV YDITNINSFK ETDKWVEHVR AERGDDVVIM LLGNKIDLDN KRQVTAEQGE EKSRNLNVMF IETSAKTGYN VKKLFRRVAS ALLSTRTSPP PKEGTVEIEL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com