

Product Data Sheet

Recombinant Human RABL6 Protein

Catalog #	Source	Reactivity	Applications
CST1320	E. coli	H	
Description	Recombinant Human RABL6 Protein is produced by E. coli expression system and the target gene encoding Ser11-Pro330 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RABL6		
Alternative Names	C9orf86; PARF; Rab-like protein 6; GTP-binding protein Parf; Partner of ARF; Rab-like protein 1; RBEL1		
Entrez Gene	55684		
SwissProt	Q3YEC7		
Purity			
Chemical Structure	SDQAPGRDKN IPAGLQSMNQ ALQRRFAKGV QYNMKIVIRG DRNTGKTALW HRLQGRPFVE EYIPTQEIQV TSIHWSYKTT DDIVKVEVWD VVDKGKCKKR GDGLKMENDP QEAESMALD AEFLDVYKNC NGVMMFDIT KQWTFNYILR ELPKVPHTVP VCVLGNYRDM GEHRVILPDD VRDFIDNLDR PPGSSYFRIA ESSMKNSFGL KYLHKFFNIP FLQLQRETL RQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com