

Product Data Sheet

Recombinant Human VMAC Protein

Catalog #	Source	Reactivity	Applications
CST1321	E. coli	H	
Description	Recombinant Human VMAC Protein is produced by E. coli expression system and the target gene encoding Leu51-Leu160 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	VMAC		
Alternative Names	Vimentin-type intermediate filament-associated coiled-coil protein		
Entrez Gene	400673		
SwissProt	Q2NL98		
Purity			
Chemical Structure	LRAALDELGR AKDREIATLQ EQLMTSEATV HSLQATVHQR DELIRQLQPR AELLQDICRR RPPLAGLLDA LAEAERLGPL PASDPGHPPP GGGPPLDNS TGEEADRDHL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com