

Product Data Sheet

Recombinant Human CIBAR1 Protein

Catalog #	Source	Reactivity	Applications
CST1335	E. coli	H	
Description	Recombinant Human CIBAR1 Protein is produced by E. coli expression system and the target gene encoding Leu151-Asp280 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FAM92A1		
Alternative Names	Protein FAM92A1		
Entrez Gene	137392		
SwissProt	A1XBS5		
Purity			
Chemical Structure	LQRAAMDASR TSRHLEETIN NFERQKMKDI KTIFSEFITI EMLFHGKALE VYTAAYQNIQ NIDEDEDLEV FRNSLYAPDY SSRLDIVRAN SKSPLQRSLS AKCVSGTGQV STCRLRKDQQ AEDDEDDELD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com