

Product Data Sheet

Recombinant Human C9orf78 Protein

Catalog #	Source	Reactivity	Applications
CST1343	E. coli	H	
Description	Recombinant Human C9orf78 Protein is produced by E. coli expression system and the target gene encoding Glu61-Lys240 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	C9orf78		
Alternative Names	HCA59; Uncharacterized protein C9orf78; Hepatocellular carcinoma-associated antigen 59		
Entrez Gene	51759		
SwissProt	Q9NZ63		
Purity			
Chemical Structure	ETTLVDDPFQ MKTGGMVDMK KLKERGKDKI SEEEDLHLGT SFS AETNRRD EDADMMKYIE TELKKRK GIV EHEEQKV KPK NAEDCLYELP ENIRVSSAKK TEEMLSNQML SG IPEVDLGI DAKIKNIIST EDAKARLLAE QQNKKK DSET SFVPTNMAVN YVQHNR FYHE ELNAPIRRNK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com