

Product Data Sheet

Recombinant Human SNX16 Protein

Catalog #	Source	Reactivity	Applications
CST1347	E. coli	H	
Description	Recombinant Human SNX16 Protein is produced by E. coli expression system and the target gene encoding Val101-Glu310 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SNX16		
Alternative Names	Sorting nexin-16		
Entrez Gene	64089		
SwissProt	P57768		
Purity			
Chemical Structure	VNWEDRPSTP TILGYEVMEE RAKFTVYKIL VKKTPEESWV VFRRYTDFSR LNDKCLKEMFP GFRLLALPPKR WFKDNYNADF LEDRQLGLQA FLQNLVAHKD IANCLAVREF LCLDDPPGPF DSLEESRAFC ETLEETNYRL QKELLEKQKE MESLKKLLSE KQLHIDTLEN RIRTLSLEPE ESLDVSETEG EQILKVESSA LEVDQDVLDE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com