

Product Data Sheet

Recombinant Human ZNF788P Protein

Catalog #	Source	Reactivity	Applications
CST1354	E. coli	H	
Description	Recombinant Human ZNF788P Protein is produced by E. coli expression system and the target gene encoding Met1-Phe82 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ZNF788		
Alternative Names	Zinc finger protein 788		
Entrez Gene	388507		
SwissProt	Q6ZQV5		
Purity			
Chemical Structure	MRNMIPQDNE NPPQQGEANQ NDSVAFEDVA VNFTPDEWAL LDPSQKNLYR EVMQETLRNL ASIEVLWKRD SLKVKVISME KF		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com