

Recombinant Human PRRC2C Protein

Catalog #	Source	Reactivity	Applications
CST1359	E. coli	H	
Description	Recombinant Human PRRC2C Protein is produced by E. coli expression system and the target gene encoding Arg471-Gln690 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	PRRC2C		
Alternative Names	BAT2D1; BAT2L2; KIAA1096; XTP2; Protein PRRC2C; BAT2 domain-containing protein 1; HBV X-transactivated gene 2 protein; HBV XAg-transactivated protein 2; HLA-B-associated transcript 2-like 2; Proline-rich and coiled-coil-containing protein 2C		
Entrez Gene	23215		
SwissProt	Q9Y520		
Purity			
Chemical Structure	RRMEEQRKAA CAEKLKRLDE KLGILEKQPS PEEIRERERE KEREREKELE KEQEQUEREKE REKDRERQQE KEKELEKEQE KQREMEKERK QEKEKELERQ KEKEKELQKM KEQEKECELE KEREKLEEKI EPREP NLEPM VEKQSENSC NKEEEPVFTR QDSNRSEKEA TPVVHETEPE SGSQPRPAVL SGYFKQFQKS LPPRFQRQQE QMKQQQWQQQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com