

Product Data Sheet

Recombinant Human KATNIP Protein

Catalog #	Source	Reactivity	Applications
CST1360	E. coli	H	
Description	Recombinant Human KATNIP Protein is produced by E. coli expression system and the target gene encoding Ala611-Thr960 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	KIAA0556		
Alternative Names	Uncharacterized protein KIAA0556		
Entrez Gene	23247		
SwissProt	O60303		
Purity			
Chemical Structure	APADHSILVD QKNEKSEQL EAMNAHSEES KGTHEMAGAS GDKELGLGCS PPAETLADAK LSSQGNVSGK RKNSTNCRKD SLSQLEEYLR LSAVPTSMGD MPSAPATSPP VKCPPVHEEP SLIQQLENLM GRKICEPPGK TPSWLQPSPT GKDRKQGGRK PKPLWLSPEK PLAWKGRGPS DDVIGEGPGE TEARDKGLRH EPGWGTSSRV NTKERPQRAT TKVHSDSDI FN		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com