

Product Data Sheet

Recombinant Human JADE3 Protein

Catalog #	Source	Reactivity	Applications
CST1361	E. coli	H	
Description	Recombinant Human JADE3 Protein is produced by E. coli expression system and the target gene encoding Pro621-Lys810 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	JADE3		
Alternative Names	KIAA0215; PHF16; Protein Jade-3; Jade family PHD finger protein 3; PHD finger protein 16		
Entrez Gene	9767		
SwissProt	Q92613		
Purity			
Chemical Structure	PAWRTPSSEC YHGQSLGKPL VLQAALHGQS SIGNGKSQPN SKFAKSNGLE GSWSGNVTQK DSSSEMFCDDQ EPVFSPHLVS QGSFRKSTVE HFSRSFKETT NRWVKNTEDL QCYVKPTKNM SPKEQFWGRQ VLRRSAGRAP YQENDGYCPD LELSDSEAES DGNKEKVRVR KDSSDRENPP HDSRRDCHGK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com