

Product Data Sheet

Recombinant Human FRYL Protein

Catalog #	Source	Reactivity	Applications
CST1362	E. coli	H	
Description	Recombinant Human FRYL Protein is produced by E. coli expression system and the target gene encoding Leu2541-Lys2860 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FRYL		
Alternative Names	AF4P12; KIAA0826; Protein furry homolog-like; ALL1-fused gene from chromosome 4p12 protein		
Entrez Gene	285527		
SwissProt	O94915		
Purity			
Chemical Structure	LQIRDETPTL EASLDNANSR LPEDTTSVLK EEHVTTFEDE GSYIIQEQQE SLVCQGILDL EETEMPEPLA PESYPESVCE EDVTLALKEL DERCEEEEEAD FSGLSQDEE EQDGFPEVQT SPLSPFLSA IIAAFQPVAY DDEEEAWRCH VNQMLSDTDG SSAVFTFHVF SRLFQTIQRK FGEITNEAVS FLGDSLQRIG TKFKSSLEVM MLCSECPTVF VDAETLMSCG LL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com