

Product Data Sheet

Recombinant Human FAM83A Protein

Catalog #	Source	Reactivity	Applications
CST1363	E. coli	H	
Description	Recombinant Human FAM83A Protein is produced by E. coli expression system and the target gene encoding Asp61-Pro310 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FAM83A		
Alternative Names	TSGP; Protein FAM83A; Tumor antigen BJ-TSA-9; Tumor-specific gene expressed in prostate protein		
Entrez Gene	84985		
SwissProt	Q86UY5		
Purity			
Chemical Structure	DFLSSVEAQY IQAQAREPPC PPDTLGGAEA GPKGLDSSSL QSGTYFPVAS EGSEPALLHS WASAEKPYLK EKSSATVYFQ TVKHNNIRDL VRRICITRTSQ VLVILMDVFT DVEIFCDILE AANKRGVFVC VLLDQGGVKL FQEMCDKVQI SDSHLKNISI RSVEGEIYCA KSGRKFAGQI REKFIISDWR FVLSGSYSFT WLCGHVHRNI LSKFTGQAVE LFDEEFRHLY AS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com