

Product Data Sheet

Recombinant Human HYCC2 Protein

Catalog #	Source	Reactivity	Applications
CST1364	E. coli	H	
Description	Recombinant Human HYCC2 Protein is produced by E. coli expression system and the target gene encoding Ser341-Pro520 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FAM126B		
Alternative Names	Protein FAM126B		
Entrez Gene	285172		
SwissProt	Q8IXS8		
Purity			
Chemical Structure	SVNLNDADEG FSSGASLSSQ PIGTKPSSSS QRGSLRKVAT GRSKDKETA SAIKSSESPR DSVVRKQYVQ QPTDLSVDSV ELTPMKKHLS LPAGQVVPKI NSLSLIRTAS ASSSKSFDYV NGSQASTSIG VGTEGGTNLA ANNANRYSTV SLQEDRLGQA GEGKELLSPG APLTKQSRSP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com