

Product Data Sheet

Recombinant Human DOP1B Protein

Catalog #	Source	Reactivity	Applications
CST1365	E. coli	H	
Description	Recombinant Human DOP1B Protein is produced by E. coli expression system and the target gene encoding Glu561-Asn680 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	DOPEY2		
Alternative Names	C21orf5; KIAA0933; Protein dopey-2		
Entrez Gene	9980		
SwissProt	Q9Y3R5		
Purity			
Chemical Structure	ENGKIILETK AVIPGDEDAS FPPLKSEDSG IGLSASSPEL SEHLRVPRVS LERDDVWKKG GSMQRTFLCI QELIANFASK NIFGVQLTAS GEESKSEPA GKRDRDGTQS LAANDSSRKN		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com