

Product Data Sheet

Recombinant Human CFAP418 Protein

Catalog #	Source	Reactivity	Applications
CST1366	E. coli	H	
Description	Recombinant Human CFAP418 Protein is produced by E. coli expression system and the target gene encoding Met1-His207 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	C8orf37		
Alternative Names	Protein C8orf37		
Entrez Gene	157657		
SwissProt	Q96NL8		
Purity			
Chemical Structure	MAEDLDELDD EVESKFCTPD LLRRGMVEQP KGCGGGTHSS DRNQAKAKET LRSTETFKKE DDLDSLNI EI LEEPNLDDKKP SKLKS KSSGN TSVRASIEGL GKSCSPVYLG GSSIPCGIGT NISWRACDHL RCIACDFLVV SYDDYMWDKS CDYLFFRNNM PEFHKLKAKL IKKKGTRAYA CQCSWRTIEE VTDLQTDHQL RWVCGKH		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com