

Product Data Sheet

Recombinant Human TRIML1 Protein

Catalog #	Source	Reactivity	Applications
CST1367	E. coli	H	
Description	Recombinant Human TRIML1 Protein is produced by E. coli expression system and the target gene encoding Met1-Val468 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TRIML1		
Alternative Names	RNF209; Probable E3 ubiquitin-protein ligase TRIML1; RING finger protein 209; Tripartite motif family-like protein 1		
Entrez Gene	339976		
SwissProt	Q8N9V2		
Purity			
Chemical Structure	MSTADLMENL REELTCFICL DYFSSPVTE CGHSFCLVCL LRSWEEHNTPLSCPECWRTL EGPHFQSNER LGRLASIARQ LRSQVLQSED EQGSYGRMPT TAKALSDDEQ GGSAFVAQSH GANRVHLSSE AEEHHREKLQ EILNLLRVRR KEAQAVLTHERKVRKLCQEE TKTCKQVVVS EYMKMHQFLK EEEQLQLQLL EQEEKENMRK LRNNEIKLTQ QIRSLSKMIA QI		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com