

Product Data Sheet

Recombinant Human PLEKHA3 Protein

Catalog #	Source	Reactivity	Applications
CST1369	E. coli	H	
Description	Recombinant Human PLEKHA3 Protein is produced by E. coli expression system and the target gene encoding Pro191-Glu290 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	PLEKHA3		
Alternative Names	FAPP1; Pleckstrin homology domain-containing family A member 3; PH domain-containing family A member 3; Phosphatidylinositol-four-phosphate adapter protein 1; FAPP-1; Phosphoinositol 4-phosphate adapter protein 1		
Entrez Gene	65977		
SwissProt	Q9HB20		
Purity			
Chemical Structure	PDPLVSPVSP SPVQMMKRSV SHPGSCSSER SSHSIKEPVS TLHRLSQRRR RTYSDTDSCS DIPLEDPDRP VHCSKNTLNG DLASATIPPEE SRLMAKKQSE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com