

Product Data Sheet

Recombinant Human OSBPL11 Protein

Catalog #	Source	Reactivity	Applications
CST1370	E. coli	H	
Description	Recombinant Human OSBPL11 Protein is produced by E. coli expression system and the target gene encoding Ser281-Glu490 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	OSBPL11		
Alternative Names	ORP11; OSBP12; Oxysterol-binding protein-related protein 11; ORP-11; OSBP-related protein 11		
Entrez Gene	114885		
SwissProt	Q9BXB4		
Purity			
Chemical Structure	SHQKGSLPSG TTIEWLEPKI SLSNHYKNGA DQPFATDQSK PVAVPEEQPV AESGLLAREP EEINADDEIE DTCDHKEDDL GAVEEQRSVI LHLLSQLKLG MDLTRVVLPT FILEKRSLL MYADFMSHPD LFIAITNGAT AEDRMIRFVE YYLTSFHEGR KGAIKKPYN PIIGETFHCS WKMPKSEVAS SVFSSSTQG VTNHAPLSGE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com