

Product Data Sheet

Recombinant Human TIMM13 Protein

Catalog #	Source	Reactivity	Applications
CST1376	E. coli	H	
Description	Recombinant Human TIMM13 Protein is produced by E. coli expression system and the target gene encoding Met1-Met95 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TIMM13		
Alternative Names	TIM13B; TIMM13A; TIMM13B; Mitochondrial import inner membrane translocase subunit Tim13		
Entrez Gene	26517		
SwissProt	Q9Y5L4		
Purity			
Chemical Structure	MEGGFGSDFG GSGSGKLDPG LIMEQVKVQI AVANAQELLQ RMTDKCFRKC IGKPGGSLDN SEQKCIAMCM DRYMDAWNTV SRAYNSRLQR ERANM		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com