

Product Data Sheet

Recombinant Human LRGUK Protein

Catalog #	Source	Reactivity	Applications
CST1382	E. coli	H	
Description	Recombinant Human LRGUK Protein is produced by E. coli expression system and the target gene encoding Glu101-Leu410 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LRGUK		
Alternative Names	Leucine-rich repeat and guanylate kinase domain-containing protein		
Entrez Gene	136332		
SwissProt	Q96M69		
Purity			
Chemical Structure	EFDGVLREEA VAKALHHLGR SGSGTEQVYL NLTLSGCNLI DVSILCGYVH LQKLDLSANK IEDLSCVSCM PYLLELNASQ NNLTTFNFK PPKNLKKADF SHNQISEICD LSAYHALTKL ILDGNEIEEI SGLEMCNNLI HLSLANNKIT TINGLNKLPI KILCLSNNQI EMITGLEDLK ALQNLDLSHN QISSLQGLN HDLLEVINLE DNKIAELREI EYIKNLPILR VL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com