

Product Data Sheet

Recombinant Human LURAP1 Protein

Catalog #	Source	Reactivity	Applications
CST1383	E. coli	H	
Description	Recombinant Human LURAP1 Protein is produced by E. coli expression system and the target gene encoding Asp11-Pro210 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LURAP1		
Alternative Names	C1orf190; LRAP35A; LRP35A; Leucine rich adaptor protein 1; Leucine repeat adapter protein 35A		
Entrez Gene	541468		
SwissProt	Q96LR2		
Purity			
Chemical Structure	DLRDVEGKVG RKTPEGLLRG LRGECELGTS GALLLPGASS TGHDLGDKIM ALKMELAYLR AIDVKILQQL VTLNEGIEAV RWLLEERGTL TSHCSSLTSS QYSLTGGSPG RSRRGSWDSL PDTSTTDRLD SVSIGSFOLDT VAPSELDEQG PPGAPRSEMD WAKVIAGGER ARTEVDVAAT RLGSLRAVWK PPGERLQGGP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com