

Product Data Sheet

Recombinant Human KIF16B Protein

Catalog #	Source	Reactivity	Applications
CST1384	E. coli	H	
Description	Recombinant Human KIF16B Protein is produced by E. coli expression system and the target gene encoding Glu691-Arg830 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	KIF16B		
Alternative Names	C20orf23; KIAA1590; SNX23; Kinesin-like protein KIF16B; Sorting nexin-23		
Entrez Gene	55614		
SwissProt	Q96L93		
Purity			
Chemical Structure	ELQKKRQEEE TFLRVQEELQ RLKELNNNEK AEKFQIFQEL DQLQKEKDEQ YAKLELEKKR LEEQEKEQVM LVAHLEEQLR EKQEMIQLLR RGEVQWVEEE KRDLEGIRES LLRVKEARAG GDEEDGEELEK AQLRFFFEFKR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com