

Recombinant Human ITSN2 Protein

Catalog #	Source	Reactivity	Applications
CST1386	E. coli	H	
Description	Recombinant Human ITSN2 Protein is produced by E. coli expression system and the target gene encoding Arg381-Glu750 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ITSN2		
Alternative Names	KIAA1256; SH3D1B; SWAP; Intersectin-2; SH3 domain-containing protein 1B; SH3P18; SH3P18-like WASP-associated protein		
Entrez Gene	50618		
SwissProt	Q9NZM3		
Purity			
Chemical Structure	RRQALMEQQQ REAERKAQKE KEEWERKQRE LQEQEWKKQL ELEKRLEKQR ELERQREEER RKDIERREAA KQELERQRR L EWERIRRRQEL LNQKNREQEE IVRLNSKKKN LHLELEALNG KHQQISGRLQ DVRLKKQTQK TELEVLDKQC DLEIMEIKQL QQELQEYQNK LIYLVPEKQL LNERIKNMQF SNTPDSGVSL LHKKSLEKEE LCQRLKEQLD ALEKETASKL SE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com