

Product Data Sheet

Recombinant Human INTS6L Protein

Catalog #	Source	Reactivity	Applications
CST1387	E. coli	H	
Description	Recombinant Human INTS6L Protein is produced by E. coli expression system and the target gene encoding Ala671-Thr820 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	DDX26B		
Alternative Names	Protein DDX26B		
Entrez Gene	203522		
SwissProt	Q5JSJ4		
Purity			
Chemical Structure	ATIIHDGHEE KMENGQITPD GFLSKSAPSE LINMTGDLMP PNQVDSLSD FTSLSKDGLI QKPGSNAFVG GAKNCSLSVD DQKDPVASTL GAMPNTLQIT PAMAQGINAD IKHQLMKEVR KFGKRYERIF ILLEEVQGPL EMKKQFVEFT		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com