

Product Data Sheet

Recombinant Human GARNL3 Protein

Catalog #	Source	Reactivity	Applications
CST1390	E. coli	H	
Description	Recombinant Human GARNL3 Protein is produced by E. coli expression system and the target gene encoding Ser811-Gly990 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GARNL3		
Alternative Names	GTPase-activating Rap/Ran-GAP domain-like protein 3		
Entrez Gene	84253		
SwissProt	Q5VWW2		
Purity			
Chemical Structure	SSGGSSKGAS ARNSPQTTPG RDTPVFPSSL GEGEIQSKNL YKIPLRNLVG RSIERPLKSP LVSKVITPPT PISVGLAAIP VTHSLSLSRM EIKEIASRTR RELLGLSDEG GPKSEGAPKA KSKPRKRLEE SQGGPKPGAV RSSSSDRIPS GSLESASTSE ANPEGHSASS DQDPVADREG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com