

Product Data Sheet

Recombinant Human GOLGA8CP Protein

Catalog #	Source	Reactivity	Applications
CST1391	E. coli	H	
Description	Recombinant Human GOLGA8CP Protein is produced by E. coli expression system and the target gene encoding His281-Val510 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GOLGA8CP		
Alternative Names	GOLGA8C; Golgin subfamily A member 8C		
Entrez Gene	729786		
SwissProt	A6NN73		
Purity			
Chemical Structure	HQMAEPLPPE PPAVPSEVEL QHLRKELERV AGALQAQVEY NQRISLLNEG QKERLREQEE RLQEQQERLR EQEERLQQLA EPQNSFKELN NENKSVLQLE QQVKELQEKL GKRLAHPVAS AQKEPEAAVP APGPGGESSG FMDHLEEKAD LSELVEKEEL GFFQYYRERC HQKVYHPITK PGGSAKDAAP GGGHHQAGPG QGGDEGEAAG AAGDGVAAGG DYKGHSKFLV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com