

Product Data Sheet

Recombinant Human EXOC3L4 Protein

Catalog #	Source	Reactivity	Applications
CST1392	E. coli	H	
Description	Recombinant Human EXOC3L4 Protein is produced by E. coli expression system and the target gene encoding Pro11-Leu340 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	EXOC3L4		
Alternative Names	C14orf73; Exocyst complex component 3-like protein 4		
Entrez Gene	91828		
SwissProt	Q17RC7		
Purity			
Chemical Structure	PELQSPKEAE EPQTPAQGSR RTSSRKEPNA HRKDGTRLGL GSLRQAFSRA SQRALTQVSK EDTGLFRRSS CSLFRSFRQA LNDGPATGHS QATPEVPSGV MNGVSQQAST GAASEELKPE AEGKSVADLI TERQLLAAFE QLLRLETLV AEKASRTFEQ DPTAFARRAM DVCLLYDGLA AEIGAIVRET LDSDGVDAAA LAELARVVSA EEEAHPSPPD DGDFLRTPRR WR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com