

Product Data Sheet

Recombinant Human EPS8L1 Protein

Catalog #	Source	Reactivity	Applications
CST1393	E. coli	H	
Description	Recombinant Human EPS8L1 Protein is produced by E. coli expression system and the target gene encoding Gly361-Met720 is expressed with N-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	EPS8L1		
Alternative Names	DRC3; EPS8R1; Epidermal growth factor receptor kinase substrate 8-like protein 1; EPS8-like protein 1; Epidermal growth factor receptor pathway substrate 8-related protein 1; EPS8-related protein 1		
Entrez Gene	54869		
SwissProt	Q8TE68		
Purity			
Chemical Structure	GGPEFASSVR RPHLTSDAVA LLRDNVTPRE NELWTSLGDS WTRPGLELSP EEGPPYRPEF FSGWEPPVTD PQSRAWEDPV EKQLQHERRR RQQSAPQVAV NGHRLPEPES EPQLESETAG KWVLCNYDFQ ARNSSELSVK QRDVLEVLDD SRKWWKVRDP AGQEGYVPYN ILTPYPGPRL HHSQSPARSL NSTPPPPPAP APAPPPALAR PRWDRPRWDS CDSLNGLDPS EK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com