

Product Data Sheet

Recombinant Human DLGAP4 Protein

Catalog #	Source	Reactivity	Applications
CST1395	E. coli	H	
Description	Recombinant Human DLGAP4 Protein is produced by E. coli expression system and the target gene encoding His511-Pro660 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	DLGAP4		
Alternative Names	DAP4; KIAA0964; SAPAP4; Disks large-associated protein 4; DAP-4; PSD-95/SAP90-binding protein 4; SAP90/PSD-95-associated protein 4; SAPAP-4		
Entrez Gene	22839		
SwissProt	Q9Y2H0		
Purity			
Chemical Structure	HSYLRAIQAG CSQEEDSVSL QSLPPPSTG SLSNSRTLPS SSCLVAYKKT PPPVPPRTTS KPFISVTVQS STESAQDTYL DSQDHKSEVT SQSGLSNSSD SLDSSTRPPS VTRGGVAPAP EAPEPPPKHA ALKSEQGTLT SESHPEAAP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com