

Product Data Sheet

Recombinant Human DELE1 Protein

Catalog #	Source	Reactivity	Applications
CST1396	E. coli	H	
Description	Recombinant Human DELE1 Protein is produced by E. coli expression system and the target gene encoding Gln141-Asp260 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	KIAA0141		
Alternative Names	DELE; Death ligand signal enhancer		
Entrez Gene	9812		
SwissProt	Q14154		
Purity			
Chemical Structure	QHILPSDGP APRHTGLREP RLGQEEASQA PRNFSHNSLR GARPQDPSEE GPGDFGFLHA SSSIESEAKP AQPQPTGEKE QDKSKTSLSE EAVTSIQQLF QLSVSI AFNF LGTENMKSGD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com