

Recombinant Human PITPNC1 Protein

Catalog #	Source	Reactivity	Applications
CST1397	E. coli	H	
Description	Recombinant Human PITPNC1 Protein is produced by E. coli expression system and the target gene encoding Thr241-Lys330 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	PITPNC1		
Alternative Names	Cytoplasmic phosphatidylinositol transfer protein 1; Mammalian rdgB homolog beta; M-rdgB beta; MrdgBbeta; Retinal degeneration B homolog beta; RdgBbeta		
Entrez Gene	26207		
SwissProt	Q9UKF7		
Purity			
Chemical Structure	TQEATNKKIG IFPPAISISS IPLLSSVRS APSSAPSTPL STDAPEFLSV PKDRPRKKSA PETLTLPDPE KKATLNLPGM HSSDKPCRPK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com