

Product Data Sheet

Recombinant Human CTDSPL2 Protein

Catalog #	Source	Reactivity	Applications
CST1399	E. coli	H	
Description	Recombinant Human CTDSPL2 Protein is produced by E. coli expression system and the target gene encoding Gln11-Glu270 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CTDSPL2		
Alternative Names	CTD small phosphatase-like protein 2; CTDSP-like 2		
Entrez Gene	51496		
SwissProt	Q05D32		
Purity			
Chemical Structure	QSNQIQTQRT ARAKRKYSEV DDSLPSGG EK PSKNETGLLS SIKKFIKGST PKEERENPSK RSRIERDIDN NLITSTPRAG EKPNKQISRV RRKSQVNGEA GSYEMTNQHV KQNGKLEDNP SSGSPPRTTL LGTIFSPVFN FFSPANKNGT SGSDSPGQAV EAEIIVKQLD MEQVDEITTS TTTSTNGAAY SNQAVQVRPS LNNGLEEAEE TVNRDIPPLT APVTPDSGYS SA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com