

Product Data Sheet

Recombinant Human CNTROB Protein

Catalog #	Source	Reactivity	Applications
CST1400	E. coli	H	
Description	Recombinant Human CNTROB Protein is produced by E. coli expression system and the target gene encoding Pro11-Glu400 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CNTROB		
Alternative Names	LIP8; Centrobin; Centrosomal BRCA2-interacting protein; LYST-interacting protein 8		
Entrez Gene	116840		
SwissProt	Q8N137		
Purity			
Chemical Structure	PLGAEDLLSD SSEPPGLNQV SSEVTSQLYA SLRLSRQAEA TARAQLYLPS TSPHEGLDG FAQELSRSLV VGLEKNLKKK DGSKHIFEME SVRGQLQTML QTSRDTAYRD PLIPGAGSER REEDSFDSDS TATLLNTRPL QDLSPSSSAQ ALEELFPRYT SLRPGPPLNP PDFQGLRDAL DSEHTRRKHC ERHIQSLQTR VLELQQQLAV AVAADRKKDT MIEQLDKTLA RV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com