

Product Data Sheet

Recombinant Human CDC20B Protein

Catalog #	Source	Reactivity	Applications
CST1401	E. coli	H	
Description	Recombinant Human CDC20B Protein is produced by E. coli expression system and the target gene encoding Asp31-Gly180 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CDC20B		
Alternative Names	Cell division cycle protein 20 homolog B		
Entrez Gene	166979		
SwissProt	Q86Y33		
Purity			
Chemical Structure	DLKQKRSQDS ANVLDSVNAT YSDFKSNFAK RLSAEVPVAS SPITTRWQQS QTRALSSDSF GEEQSTTYLP EASGSVLKTP PEKETLTLGS RKEQLKTPSK GISETSNSAL HFCKAPHAMD RDWKESVASK GQKCLKQLFV TQNVVQQANG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com