

Product Data Sheet

Recombinant Human CAMSAP3 Protein

Catalog #	Source	Reactivity	Applications
CST1402	E. coli	H	
Description	Recombinant Human CAMSAP3 Protein is produced by E. coli expression system and the target gene encoding Asp861-Arg1090 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CAMSAP3		
Alternative Names	KIAA1543; Calmodulin-regulated spectrin-associated protein 3; Protein Nezha		
Entrez Gene	57662		
SwissProt	Q9P1Y5		
Purity			
Chemical Structure	DGSPAGAEDS LEEESASSEGE PRVGLGFFYK DEDKPEDEMA QKRASLLERQ QRRAEERARR KQWQEVEKEQ RREEAARLAQ EEAPGPAPLV SAVPMATPAP AARAPAEVEV GPRKGDFTTRQ EYERRAQLKL MDDLDKVLRP RAAGSGGPGR GGRRATRPRS GCCDDSALAR SPARGLLGSR LSKIYSQSTL SLSTVANEAH NNLGVKRPTS RAPSPSGLMS PSRLPGSRER		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com