

Product Data Sheet

Recombinant Human BORCS6 Protein

Catalog #	Source	Reactivity	Applications
CST1403	E. coli	H	
Description	Recombinant Human BORCS6 Protein is produced by E. coli expression system and the target gene encoding Glu11-Ala230 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	C17orf59		
Alternative Names	Uncharacterized protein C17orf59		
Entrez Gene	54785		
SwissProt	Q96GS4		
Purity			
Chemical Structure	ETDLLAVAEH QALVFGGGPG RTSSEPPAGL RVSGEEETEN VGGANRHPRT SPKTSSCGVV HRPEREALEN EPGPQGTLTG AGSRRGAPGA EHEPSLSSRH KNPAPPEGKP SSGRDCRRGG PGGGMDVEQQ EEEDNDEEAA AGSRAGRSFS SRLQDSRSLD GLSEACGGAG SSGSAESGAG GRRRATISSP LELEGTVSRH GDLTHFVANN LQLKIRLSGA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com