

Product Data Sheet

Recombinant Human ANAPC7 Protein

Catalog #	Source	Reactivity	Applications
CST1405	E. coli	H	
Description	Recombinant Human ANAPC7 Protein is produced by E. coli expression system and the target gene encoding Glu391-Gly540 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ANAPC7		
Alternative Names	APC7; Anaphase-promoting complex subunit 7; APC7; Cyclosome subunit 7		
Entrez Gene	51434		
SwissProt	Q9UJX3		
Purity			
Chemical Structure	EAMVMANNVY KTLGANAQTL TLLATVCLED PVTQEAKAKTL LDKALTQRPD YIKAVVKKAE LLSREQKYED GIALLRNALA NQSDCVLHRI LGDFLVAVNE YQEAMDQYSI ALSLPNDQK SLEGMQKMEK EESPTDATQE EDVDDMEGSG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com