

Product Data Sheet

Recombinant Human ASPG Protein

Catalog #	Source	Reactivity	Applications
CST1406	E. coli	H	
Description	Recombinant Human ASPG Protein is produced by E. coli expression system and the target gene encoding Thr361-Leu570 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ASPG		
Alternative Names	C14orf76; 60 kDa lysophospholipase		
Entrez Gene	374569		
SwissProt	Q86U10		
Purity			
Chemical Structure	TPPSVEERRP SLQGNTLGGG VSWLLSLSGS QEADALRNAL VPSLACAAAH AGDVEALQAL VELGSDLGLV DFNGQTPLHA AARGGHTEAV TMLLQRGVDV NTRDTDGFSP LLLAVRGRHP GVIGLLREAG ASLSTQEELE AGTELCRLAY RADLEGLQVW WQAGADLGQP GYDGHSALHV AEAAGNLAVV AFLQSLEGAV GAQAPCPEVL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com