

Product Data Sheet

Recombinant Human C10orf90 Protein

Catalog #	Source	Reactivity	Applications
CST1408	E. coli	H	
Description	Recombinant Human C10orf90 Protein is produced by E. coli expression system and the target gene encoding Val431-Pro660 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	C10orf90		
Alternative Names	FATS; Centrosomal protein C10orf90; Fragile-site associated tumor suppressor homolog		
Entrez Gene	118611		
SwissProt	Q96M02		
Purity			
Chemical Structure	VCMTVSAPPV EQKPTRHFLP IGDSSPSDDC LSRDLSEPTE RRHQSFLLKPR ILFPGFLCPL QDVCSALQED NGVQIESKFP KGDYTCDDL VVKIKECKKSE DPTTPEPSPA APSPAPRDGA GSPGLSEDCS ESQQTPARSL TLQEALVRK PQFISRSQER LKKLEHMMVQQ RKAQRKEDLR QKQSLPIRT SKKQFTIPHP LSDNLFKPKE RCISEKEMHM RSKRIYDNL P		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com