

Product Data Sheet

Recombinant Human UBE4A Protein

Catalog #	Source	Reactivity	Applications
CST1409	E. coli	H	
Description	Recombinant Human UBE4A Protein is produced by E. coli expression system and the target gene encoding Asp41-Ala440 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	UBE4A		
Alternative Names	KIAA0126; Ubiquitin conjugation factor E4 A		
Entrez Gene	9354		
SwissProt	Q14139		
Purity			
Chemical Structure	DELPAASPDDSDNSVSESLDEFDYSVAEISR SFRSQQEICEQLNINHMIQR IFLITLDNSD PSLKSGNGIP SRCVYLEEMA VELEDQDQWLD MSNVEQALFA RLLQDPGNH LINMTSSTTL NLSADRDAGE RHIFCYLYSC FQRAKEEITK VPENLLPFAV QCRNLTVSNT RTVLLTPEIY VDQNIHEQLV DLMLEAIQGA HFEDVTEFLE EVIEALILDE EVRTFPEVMI PV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com