

Product Data Sheet

Recombinant Human TUT7 Protein

Catalog #	Source	Reactivity	Applications
CST1411	E. coli	H	
Description	Recombinant Human TUT7 Protein is produced by E. coli expression system and the target gene encoding His661-Lys1010 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ZCCHC6		
Alternative Names	HS2; KIAA1711; TUT7; Terminal uridylyltransferase 7; TUTase 7; Zinc finger CCHC domain-containing protein 6		
Entrez Gene	79670		
SwissProt	Q5VYS8		
Purity			
Chemical Structure	HPDVQTKDDK LKNSVLAQGP GATSSAANTC KVQPLTLKET AESFGSPPE EMGNEHISVH PENSDCIQAD VNSDDYKGDV VYHPETGRKN EKEKVGKRGK HLLTVDQKRG EHVVCSTRN NESESTLDLE GFQNPTAKEC EGLATLDNKA DLDGESTEGT EELEDNLNHF THSVQGQTSE MIPSEDEEED DEEEEEEEEP RLTIHQREDE DGMANEDELD NTYTGSQDED AL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com