

Product Data Sheet

Recombinant Human ARHGAP29 Protein

Catalog #	Source	Reactivity	Applications
CST1412	E. coli	H	
Description	Recombinant Human ARHGAP29 Protein is produced by E. coli expression system and the target gene encoding Leu911-Phe1250 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ARHGAP29		
Alternative Names	PARG1; Rho GTPase-activating protein 29; PTPL1-associated RhoGAP protein 1; Rho-type GTPase-activating protein 29		
Entrez Gene	9411		
SwissProt	Q52LW3		
Purity			
Chemical Structure	LLSPEERDIE RSMKSLFFSS KEDIHTSESE SKIFERATSF EESERKQNAL GKCDACLSDK AQLLLDQAE SASQKIEDGK TPKPLSLKSD RSTNNVERHT PRTKIRPVSL PVDRLLLASP PNERNGRNMG NVNLDKFCKN PAFEGVNRKD AATTVCSEFN GFDQQTQKI QDKQYEQNSL TAKTTMIMPS ALQEKGVTTT LQISGDHSIN ATQPSKPYAE PVRSVREASE RR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com