

Product Data Sheet

Recombinant Human SEC16A Protein

Catalog #	Source	Reactivity	Applications
CST1414	E. coli	H	
Description	Recombinant Human SEC16A Protein is produced by E. coli expression system and the target gene encoding Gly501-Ala650 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SEC16A		
Alternative Names	KIAA0310; SEC16; SEC16L; Protein transport protein Sec16A; SEC16 homolog A		
Entrez Gene	9919		
SwissProt	O15027		
Purity			
Chemical Structure	GAVCHTGAPD ATLHTVHPDS VSSSYSSRSRSH GRLSGSARPQ ELVGTFIIQQE VGKPEDEASG SFFKQIDSSP VGGETDETTV SQNYRGSVSQ PSTPSPPKPT GIFQTSANSS FEPVKSHLVG VKPFEADLAN VVGEVRETCV RQKQCRPAAA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com