

Product Data Sheet

Recombinant Human SBNO1 Protein

Catalog #	Source	Reactivity	Applications
CST1415	E. coli	H	
Description	Recombinant Human SBNO1 Protein is produced by E. coli expression system and the target gene encoding Lys751-Tyr1100 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SBNO1		
Alternative Names	MOP3; Protein strawberry notch homolog 1; Monocyte protein 3; MOP-3		
Entrez Gene	55206		
SwissProt	A3KN83		
Purity			
Chemical Structure	KNMSSGDDDD FNPFLDESNE DDENDPWLIR KDHKKNKEKK KKKSIDPDSI QSALLASGLG SKRPSFSSTP VISPAPNSTP ANSNTNSNSS LITSQDAVER AQQMKKDLLD KLEKLAEDLP PNTLDELIDE LGGPENVAEM TGRKGRVVSND DGSIYESR SELDVPVEIL NITEKQRFMD GDKNIAISE AASSGISLQA DRRAKNQRRR VHMTLELPWS ADRAIQQFGR TH		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com