

Product Data Sheet

Recombinant Human RNF123 Protein

Catalog #	Source	Reactivity	Applications
CST1423	E. coli	H	
Description	Recombinant Human RNF123 Protein is produced by E. coli expression system and the target gene encoding Glu451-Gly520 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RNF123		
Alternative Names	KPC1; E3 ubiquitin-protein ligase RNF123; Kip1 ubiquitination-promoting complex protein 1; RING finger protein 123		
Entrez Gene	63891		
SwissProt	Q5XPI4		
Purity			
Chemical Structure	ELIPTTWWPH CSSREGKEST EMKEETAER LRRRAYERGC QRLRKRIEVV EELQVQILKL LLDNKDDNGG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com